


# Mid-Range SEATBACK BAG

Your Seatback Bag has straps that allow you to attach the bag to your seatback in a variety of ways. If your vehicle has a cushioned seatback, loosen the horizontal strap and slide over the seatback, with the 2 vertical straps going over the top of the seatback. Then tighten the strap to secure the bag.




You can also hang your bag from push-handles. Slide the horizontal strap out of the buckle and back out of the vertical straps. The vertical straps can then slide over your push-handles. You can then wrap the horizontal strap around your seatback to keep the back firmly in place, but this is not necessary. If you choose not to use the horizontal strap, simply slide it back through the buckle or cut it off with scissors.

NOTE: If the bag is hanging under your push-handles lower than you would like, use the 2 plastic slide (included inside bag, use 1 for each side) to shorten the straps (SEE BELOW.)



Fold strap over, then insert through plastic slide (left), creating a new "loop", then hook over push-handles (right)



***We strongly encourage you to call us at (800) 795-2392 if you have any problems, questions, or comments. We are always in development of new products for wheelchair, scooter, crutch, and walker users...so let us know if there's something special you need!***